



## Run Info

Host Name	GXB03020 (localhost)
Experiment Name	ReadUntil_38Kbp_LowtoHigh_MmEnr_15042021
Sample ID	ReadUntil_38Kbp_LowtoHigh_MmEnr_15042021
Run ID	5a06b53c-cb93-4051-a39a-d7a11d821f83
Flow Cell Id	FAP14753
Start Time	April 15, 18:38
Run Length	3d 0h 3m

## Run Summary

Reads Generated	1.46 M
Passed Bases	8.39 Gb
Failed Bases	607.12 Mb
Estimated Bases	9.09 Gb


## Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
Read Until	reference_files=[ ["/data/M_morganii_ref.fasta"],filter_type=enrich,first_channel=1,last_channel=256
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=[ ["/data/the7references.fasta"]
Read Filtering	min_qscore=7


## Versions

MinKNOW	21.02.5
MinKNOW Core	4.2.5
Bream	6.1.10
Guppy	4.3.4

Cumulative Output Reads




Cumulative Output Bases




## Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 29.22 K




## Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 28.81 K




### Read Length Histogram Estimated Bases

Estimated N50: 29.22 K




## Read Length Histogram Basecalled Bases


Estimated N50: 28.81 K




Alignment Target Hits (reads)




Alignment Target Hits (bases)




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised



Temperature History.




Bias Voltage History.







Translocation Speed



QScore



Disk Write Performance



## **Run Debug Messages**

- The sequencing run has finished, but basecalling may continue April 18, 18:41
- Mux scan for flow cell FAP14753 has found a total of 76 pores. 73 pores available for immediate sequencing April 18, 17:57
- Performing Mux Scan April 18, 17:54
- Mux scan for flow cell FAP14753 has found a total of 69 pores. 66 pores available for immediate sequencing April 18, 16:24
- Performing Mux Scan April 18, 16:22
- Mux scan for flow cell FAP14753 has found a total of 71 pores. 66 pores available for immediate sequencing April 18, 14:52
- Performing Mux Scan April 18, 14:49
- Mux scan for flow cell FAP14753 has found a total of 100 pores. 93 pores available for immediate sequencing April 18, 13:19
- Performing Mux Scan April 18, 13:17
- Mux scan for flow cell FAP14753 has found a total of 62 pores. 61 pores available for immediate sequencing April 18, 11:47
- Performing Mux Scan April 18, 11:45
- Mux scan for flow cell FAP14753 has found a total of 71 pores. 69 pores available for immediate sequencing April 18, 10:14
- Performing Mux Scan April 18, 10:12
- Mux scan for flow cell FAP14753 has found a total of 63 pores. 62 pores available for immediate sequencing April 18, 08:42
- Performing Mux Scan April 18, 08:40
- Mux scan for flow cell FAP14753 has found a total of 46 pores. 45 pores available for immediate sequencing April 18, 07:09
- Performing Mux Scan April 18, 07:07
- Mux scan for flow cell FAP14753 has found a total of 70 pores. 69 pores available for immediate sequencing April 18, 05:37
- Performing Mux Scan April 18, 05:35
- Mux scan for flow cell FAP14753 has found a total of 54 pores. 52 pores available for immediate sequencing April 18, 04:05
- Performing Mux Scan April 18, 04:02
- Mux scan for flow cell FAP14753 has found a total of 53 pores. 51 pores available for immediate sequencing April 18, 02:32
- Performing Mux Scan April 18, 02:30
- Mux scan for flow cell FAP14753 has found a total of 52 pores. 52 pores available for immediate sequencing April 18, 01:00
- Performing Mux Scan April 18, 00:57
- Mux scan for flow cell FAP14753 has found a total of 61 pores. 55 pores available for immediate sequencing April 17, 23:27
- Performing Mux Scan April 17, 23:25
- Mux scan for flow cell FAP14753 has found a total of 48 pores. 47 pores available for immediate sequencing April 17, 21:55
- Performing Mux Scan April 17, 21:52
- Mux scan for flow cell FAP14753 has found a total of 77 pores. 71 pores available for immediate sequencing April 17, 20:22
- Performing Mux Scan April 17, 20:20
- Mux scan for flow cell FAP14753 has found a total of 80 pores. 74 pores available for immediate sequencing April 17, 18:50
- Performing Mux Scan April 17, 18:47

- Mux scan for flow cell FAP14753 has found a total of 74 pores. 69 pores available for immediate sequencing April 17, 17:17
- Performing Mux Scan April 17, 17:15
- Mux scan for flow cell FAP14753 has found a total of 98 pores. 86 pores available for immediate sequencing April 17, 15:45
- Performing Mux Scan April 17, 15:42
- Mux scan for flow cell FAP14753 has found a total of 79 pores. 74 pores available for immediate sequencing April 17, 14:12
- Performing Mux Scan April 17, 14:10
- Mux scan for flow cell FAP14753 has found a total of 100 pores. 94 pores available for immediate sequencing April 17, 12:40
- Performing Mux Scan April 17, 12:37
- Mux scan for flow cell FAP14753 has found a total of 108 pores. 101 pores available for immediate sequencing April 17, 11:07
- Performing Mux Scan April 17, 11:05
- Mux scan for flow cell FAP14753 has found a total of 105 pores. 96 pores available for immediate sequencing April 17, 09:35
- Performing Mux Scan April 17, 09:32
- Mux scan for flow cell FAP14753 has found a total of 94 pores. 90 pores available for immediate sequencing April 17, 08:02
- Performing Mux Scan April 17, 08:00
- Mux scan for flow cell FAP14753 has found a total of 74 pores. 71 pores available for immediate sequencing April 17, 06:30
- Performing Mux Scan April 17, 06:27
- Mux scan for flow cell FAP14753 has found a total of 123 pores. 110 pores available for immediate sequencing April 17, 04:57
- Performing Mux Scan April 17, 04:55
- Mux scan for flow cell FAP14753 has found a total of 119 pores. 102 pores available for immediate sequencing April 17, 03:25
- Performing Mux Scan April 17, 03:22
- Mux scan for flow cell FAP14753 has found a total of 124 pores. 107 pores available for immediate sequencing April 17, 01:52
- Performing Mux Scan April 17, 01:50
- Mux scan for flow cell FAP14753 has found a total of 169 pores. 131 pores available for immediate sequencing April 17, 00:19
- Performing Mux Scan April 17, 00:17
- Mux scan for flow cell FAP14753 has found a total of 322 pores. 236 pores available for immediate sequencing April 16, 22:46
- Performing Mux Scan April 16, 22:44
- Mux scan for flow cell FAP14753 has found a total of 383 pores. 269 pores available for immediate sequencing April 16, 21:13
- Performing Mux Scan April 16, 21:10
- Mux scan for flow cell FAP14753 has found a total of 376 pores. 234 pores available for immediate sequencing April 16, 19:39
- Performing Mux Scan April 16, 19:37
- Mux scan for flow cell FAP14753 has found a total of 387 pores. 229 pores available for immediate sequencing April 16, 18:06
- Performing Mux Scan April 16, 18:03
- Mux scan for flow cell FAP14753 has found a total of 476 pores. 275 pores available for immediate sequencing April 16, 16:32
- Performing Mux Scan April 16, 16:30
- Mux scan for flow cell FAP14753 has found a total of 566 pores. 323 pores available for

- immediate sequencing April 16, 14:59
- Performing Mux Scan April 16, 14:56
- Mux scan for flow cell FAP14753 has found a total of 583 pores. 316 pores available for immediate sequencing April 16, 13:25
- Performing Mux Scan April 16, 13:23
- Mux scan for flow cell FAP14753 has found a total of 578 pores. 296 pores available for immediate sequencing April 16, 11:52
- Performing Mux Scan April 16, 11:50
- Mux scan for flow cell FAP14753 has found a total of 618 pores. 309 pores available for immediate sequencing April 16, 10:18
- Performing Mux Scan April 16, 10:16
- Mux scan for flow cell FAP14753 has found a total of 680 pores. 359 pores available for immediate sequencing April 16, 08:45
- Performing Mux Scan April 16, 08:43
- Mux scan for flow cell FAP14753 has found a total of 697 pores. 348 pores available for immediate sequencing April 16, 07:11
- Performing Mux Scan April 16, 07:09
- Mux scan for flow cell FAP14753 has found a total of 704 pores. 345 pores available for immediate sequencing April 16, 05:38
- Performing Mux Scan April 16, 05:36
- Mux scan for flow cell FAP14753 has found a total of 738 pores. 367 pores available for immediate sequencing April 16, 04:05
- Performing Mux Scan April 16, 04:02
- Mux scan for flow cell FAP14753 has found a total of 793 pores. 394 pores available for immediate sequencing April 16, 02:31
- Performing Mux Scan April 16, 02:29
- Mux scan for flow cell FAP14753 has found a total of 807 pores. 395 pores available for immediate sequencing April 16, 00:58
- Performing Mux Scan April 16, 00:55
- Mux scan for flow cell FAP14753 has found a total of 854 pores. 414 pores available for immediate sequencing April 15, 23:24
- Performing Mux Scan April 15, 23:22
- Mux scan for flow cell FAP14753 has found a total of 872 pores. 424 pores available for immediate sequencing April 15, 21:51
- Performing Mux Scan April 15, 21:48
- Mux scan for flow cell FAP14753 has found a total of 889 pores. 433 pores available for immediate sequencing April 15, 20:17
- Performing Mux Scan April 15, 20:15
- Mux scan for flow cell FAP14753 has found a total of 914 pores. 440 pores available for immediate sequencing April 15, 18:44
- Performing Mux Scan April 15, 18:41
- Starting sequencing procedure April 15, 18:41
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C April 15, 18:38