

Run Info

Host Name GXB01190 (localhost)

Experiment Name ReadUntilMock_15kbSE_MwEnrich_08122020
Sample ID ReadUntilMock_15kbSE_MwEnrich_08122020

Run ID 8fd337af-bae6-48c4-8bd6-0f7ad84e9ffe

Flow Cell Id FAO52906

Start Time December 8, 17:48

Run Length 1h 13m

Run Summary

Reads Generated226.2 KPassed Bases159.57 MbFailed Bases9.94 MbEstimated Bases176.83 Mb

Run Parameters

Flow Cell Type FLO-MIN106 Kit SQK-LSK109 -185 mV Initial Bias Voltage FAST5 Output **Enabled FASTQ Output Enabled BAM Output Enabled** Active Channel Selection **Enabled** Basecalling on Specified Run Length 72 hours

reference_files=

Read Until ["/data/references/M_wisconsensis_ref.fasta"],filter_type=enrich,first

_channel=1,last_channel=512

FAST5 Reads per File 4000

FAST5 Output Options **zlib_compress,fastq,raw**

FASTQ Reads per File 4000

Mux Scan Period 1 hour 30 minutes

Reserved Pores 0 %

Basecall Model High-accuracy basecalling

Alignment reference_files=["/data/references/the7references.fasta"]

Read Filtering min_qscore=7

Versions


 MinKNOW
 20.10.6

 MinKNOW Core
 4.1.2


 Bream
 6.1.4

 Guppy
 4.2.3

Cumulative Output Reads




Cumulative Output Bases





Read Length Histogram Estimated Bases - Outliers Discarded


Estimated N50: 558




Read Length Histogram Basecalled Bases - Outliers Discarded









Read Length Histogram Basecalled Bases






Alignment Target Hits (reads)




Alignment Target

Alignment Target Hits (bases)




Alignment Target

Duty Time Grouped





Duty time Categorised






Mux Scan Grouped








Inactive


Active

Single Pore Reserved Pore Unavailable Multiple Saturated Zero Other


Temperature History




<u>Bias Voltage History</u>




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAO52906 has found a total of 965 pores. 458 pores available for immediate sequencing December 8, 17:54
- Performing Mux Scan December 8, 17:52
- Starting sequencing procedure December 8, 17:52
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C December 8, 17:48