



Run Info

Host Name	GXB01190 (localhost)
Experiment Name	ReadUntilMock_15kbSE_MmEnrich_08122020
Sample ID	ReadUntilMock_15kbSE_MmEnrich_08122020
Run ID	a3a2f40c-83a0-413b-ab2f-fb5a4f0e1a4d
Flow Cell Id	FAO52906
Start Time	December 8, 12:54
Run Length	1h 15m

Run Summary

Reads Generated	138.06 K
Passed Bases	348.4 Mb
Failed Bases	8.66 Mb
Estimated Bases	363.46 Mb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
Read Until	reference_files=["/data/references/M_morganii_ref.fasta"],filter_type=enrich,first_channel=1,last_channel=512
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=["/data/references/the7references.fasta"]
Read Filtering	min_qscore=7


Versions

MinKNOW	20.10.6
MinKNOW Core	4.1.2
Bream	6.1.4
Guppy	4.2.3

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 19.95 K




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 19.66 K




Read Length Histogram Estimated Bases


Estimated N50: 19.95 K




Alignment Target Hits (reads)




Alignment Target Hits (bases)




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History.




Bias Voltage History.




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAO52906 has found a total of 1105 pores. 478 pores available for immediate sequencing December 8, 13:00
- Performing Mux Scan December 8, 12:58
- Starting sequencing procedure December 8, 12:58
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C December 8, 12:54