

Run Info

Host Name GXb03422 (localhost)

Position X4

Experiment Name ReadUntil_ZGM_Enrich_Amuciniphila_230920
Sample ID ReadUntil_ZGM_Enrich_Amuciniphila_230920

Run ID **014e4f39-967e-480a-bb6d-828dcaecc48f**

Acquisition ID(s) c5b198cc5ac5cba429272f5fee3ddd3c7cdf306a,

e6f4367731bc587697ac48b5d7ca5c52e8e5bc8b

Flow Cell Id FAR13419

Start Time September 23, 13:42

Run Length 1h 2m

Run Summary

Reads Generated299.24 kPassed Bases412.11 MbFailed Bases28.45 MbEstimated Bases457.58 Mb

Run Parameters

FLO-MIN106 Flow Cell Type Kit SQK-LSK109 -180 mV Initial bias voltage FAST5 output **Enabled** Enabled FASTQ output Enabled BAM output Disabled Bulk file output Active channel selection **Enabled** Enabled Basecalling 72 hours Specified run length

reference_files=

Adaptive sampling ["/data/references/zymo_gut_mock/Akkermansia_muciniphila.fasta"],f

ilter_type=enrich,first_channel=1,last_channel=256

FAST5 reads per file 4000

FAST5 output options vbz_compress,fastq,raw

FASTQ reads per file 4000
FASTQ output options compress

Mux scan period 1 hour 30 minutes

Reserved pores 0 %

Basecall model High-accuracy basecalling

Alignment reference_files=

["/data/references/zymo_gut_mock/zymo_minus_two.fasta"]

Read filtering min_qscore=9


Versions

 MinKNOW
 21.05.20


 MinKNOW Core
 4.3.11

 Bream
 6.2.6

 Guppy
 5.0.13




Cumulative Output Bases



Read Length Histogram Estimated Bases - Outliers Discarded


Estimated N50: 5.58 kb


967e-480abb6d-828dcaecc48f



Read Length Histogram Basecalled Bases - Outliers Discarded


Estimated N50: 5.51 kb






Read Length Histogram Basecalled Bases


Estimated N50: 5.51 kb




Alignment Target Hits (reads)





Alignment Target Hits (bases)






Duty time Categorised












Mux Scan Categorised






Bias Voltage History











QScore



Disk Write Performance

967e-480abb6d-828dcaecc48f



ReadUntil_ZGM_Enrich_Amuciniphila_230920 ReadUntil_ZGM_Enrich_Amuciniphila_230920 014e4f39-

Run Debug Messages

967e-480abb6dFAR13419 11

828dcaecc48f

 Mux scan for flow cell FAR13419 has found a total of 1367 pores. 499 pores available for immediate sequencing September 23, 13:47

- Performing Mux Scan September 23, 13:45
- Starting sequencing procedure September 23, 13:45
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C September 23, 13:42