



Run Info

Host Name	GXb03422 (localhost)
Position	X2
Experiment Name	ReadUntil_ZGM_EcoliEnrich_B3008_23092021
Sample ID	ReadUntil_ZGM_EcoliEnrich_B3008_23092021
Run ID	21702a70-ed26-4d24-8473-a5eb92b1e81d
Acquisition ID(s)	9da1bb9b29669d3d50de0da732b688907e141af0, c5174fb226efa52d87bf22b8b152025c6ce88739
Flow Cell Id	FAR13458
Start Time	September 23, 13:05
Run Length	1h 3m

Run Summary

Reads Generated	251.52 k
Passed Bases	472.5 Mb
Failed Bases	28.56 Mb
Estimated Bases	522.29 Mb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial bias voltage	-180 mV
FAST5 output	Enabled
FASTQ output	Enabled
BAM output	Enabled
Bulk file output	Disabled
Active channel selection	Enabled
Basecalling	Enabled
Specified run length	72 hours
Adaptive sampling	reference_files=["/data/references/zymo_gut_mock/Escherichia_coli_B3008.fasta"],filter_type=enrich,first_channel=1,last_channel=256
FAST5 reads per file	4000
FAST5 output options	vbz_compress,fastq,raw
FASTQ reads per file	4000
FASTQ output options	compress
Mux scan period	1 hour 30 minutes
Reserved pores	0 %
Basecall model	High-accuracy basecalling
Alignment	reference_files=["/data/references/zymo_gut_mock/zymo_minus_two.fasta"]
Read filtering	min_qscore=9


Versions

MinKNOW	21.05.20
MinKNOW Core	4.3.11
Bream	6.2.6
Guppy	5.0.13

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 11.54 kb




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 11.24 kb




Read Length Histogram Estimated Bases

Estimated N50: 11.54 kb




Read Length Histogram Basecalled Bases


Estimated N50: 11.24 kb




Alignment Target Hits (reads)




Alignment Target Hits (bases)







Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAR13458 has found a total of 1396 pores. 500 pores available for immediate sequencing September 23, 13:11
- Performing Mux Scan September 23, 13:09
- Starting sequencing procedure September 23, 13:09
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C September 23, 13:05