



## Run Info

Host Name	GXB01190 (localhost)
Experiment Name	ReadUntil_38kb_MmDepletion_15042021
Sample ID	ReadUntil_38kb_MmDepletion_15042021
Run ID	637aed78-b222-4a41-b25f-8375a71d5c14
Flow Cell Id	FAP14879
Start Time	April 15, 10:53
Run Length	3d 0h 3m

## Run Summary

Reads Generated	1.9 M
Passed Bases	15.85 Gb
Failed Bases	1.49 Gb
Estimated Bases	17.51 Gb


## Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
Read Until	reference_files=["/data/references/M_morganii_ref.fasta"],filter_type=deplete,first_channel=1,last_channel=256
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=["/data/references/the7references.fasta"]
Read Filtering	min_qscore=7


## Versions

MinKNOW	21.02.5
MinKNOW Core	4.2.5
Bream	6.1.10
Guppy	4.3.4

**Cumulative Output Reads**




**Cumulative Output Bases**




**Read Length Histogram Estimated Bases - Outliers Discarded**

Estimated N50: 35.89 K




**Read Length Histogram Basecalled Bases - Outliers Discarded**

Estimated N50: 35.46 K




**Read Length Histogram Estimated Bases**

Estimated N50: 35.89 K




**Read Length Histogram Basecalled Bases**


Estimated N50: 35.46 K




**Alignment Target Hits (reads)**




**Alignment Target Hits (bases)**




**Duty Time Grouped**




**Duty time Categorised**




**Mux Scan Grouped**




**Mux Scan Categorised**



**Temperature History.**




**Bias Voltage History.**







**Translocation Speed**



**QScore**



**Disk Write Performance**



**Run Debug Messages**

- The sequencing run has finished, but basecalling may continue April 18, 10:56
- Mux scan for flow cell FAP14879 has found a total of 57 pores. 55 pores available for immediate sequencing April 18, 10:21
- Performing Mux Scan April 18, 10:19
- Mux scan for flow cell FAP14879 has found a total of 70 pores. 68 pores available for immediate sequencing April 18, 08:49
- Performing Mux Scan April 18, 08:46
- Mux scan for flow cell FAP14879 has found a total of 61 pores. 58 pores available for immediate sequencing April 18, 07:16
- Performing Mux Scan April 18, 07:14
- Mux scan for flow cell FAP14879 has found a total of 65 pores. 63 pores available for immediate sequencing April 18, 05:44
- Performing Mux Scan April 18, 05:41
- Mux scan for flow cell FAP14879 has found a total of 61 pores. 56 pores available for immediate sequencing April 18, 04:11
- Performing Mux Scan April 18, 04:09
- Mux scan for flow cell FAP14879 has found a total of 67 pores. 63 pores available for immediate sequencing April 18, 02:39
- Performing Mux Scan April 18, 02:36
- Mux scan for flow cell FAP14879 has found a total of 63 pores. 61 pores available for immediate sequencing April 18, 01:06
- Performing Mux Scan April 18, 01:04
- Mux scan for flow cell FAP14879 has found a total of 75 pores. 70 pores available for immediate sequencing April 17, 23:34
- Performing Mux Scan April 17, 23:32
- Mux scan for flow cell FAP14879 has found a total of 67 pores. 65 pores available for immediate sequencing April 17, 22:02
- Performing Mux Scan April 17, 21:59
- Mux scan for flow cell FAP14879 has found a total of 78 pores. 73 pores available for immediate sequencing April 17, 20:29
- Performing Mux Scan April 17, 20:27
- Mux scan for flow cell FAP14879 has found a total of 101 pores. 97 pores available for immediate sequencing April 17, 18:57
- Performing Mux Scan April 17, 18:54
- Mux scan for flow cell FAP14879 has found a total of 78 pores. 71 pores available for immediate sequencing April 17, 17:24
- Performing Mux Scan April 17, 17:22
- Mux scan for flow cell FAP14879 has found a total of 93 pores. 90 pores available for immediate sequencing April 17, 15:52
- Performing Mux Scan April 17, 15:49
- Mux scan for flow cell FAP14879 has found a total of 85 pores. 79 pores available for immediate sequencing April 17, 14:19
- Performing Mux Scan April 17, 14:17
- Mux scan for flow cell FAP14879 has found a total of 103 pores. 99 pores available for immediate sequencing April 17, 12:47
- Performing Mux Scan April 17, 12:44
- Mux scan for flow cell FAP14879 has found a total of 100 pores. 92 pores available for immediate sequencing April 17, 11:14
- Performing Mux Scan April 17, 11:12

- Mux scan for flow cell FAP14879 has found a total of 130 pores. 118 pores available for immediate sequencing April 17, 09:42
- Performing Mux Scan April 17, 09:40
- Mux scan for flow cell FAP14879 has found a total of 160 pores. 131 pores available for immediate sequencing April 17, 08:09
- Performing Mux Scan April 17, 08:06
- Mux scan for flow cell FAP14879 has found a total of 251 pores. 202 pores available for immediate sequencing April 17, 06:36
- Performing Mux Scan April 17, 06:33
- Mux scan for flow cell FAP14879 has found a total of 328 pores. 246 pores available for immediate sequencing April 17, 05:02
- Performing Mux Scan April 17, 05:00
- Mux scan for flow cell FAP14879 has found a total of 381 pores. 271 pores available for immediate sequencing April 17, 03:29
- Performing Mux Scan April 17, 03:26
- Mux scan for flow cell FAP14879 has found a total of 472 pores. 292 pores available for immediate sequencing April 17, 01:55
- Performing Mux Scan April 17, 01:53
- Mux scan for flow cell FAP14879 has found a total of 581 pores. 348 pores available for immediate sequencing April 17, 00:22
- Performing Mux Scan April 17, 00:19
- Mux scan for flow cell FAP14879 has found a total of 596 pores. 334 pores available for immediate sequencing April 16, 22:48
- Performing Mux Scan April 16, 22:46
- Mux scan for flow cell FAP14879 has found a total of 592 pores. 306 pores available for immediate sequencing April 16, 21:15
- Performing Mux Scan April 16, 21:12
- Mux scan for flow cell FAP14879 has found a total of 691 pores. 372 pores available for immediate sequencing April 16, 19:41
- Performing Mux Scan April 16, 19:39
- Mux scan for flow cell FAP14879 has found a total of 658 pores. 326 pores available for immediate sequencing April 16, 18:08
- Performing Mux Scan April 16, 18:05
- Mux scan for flow cell FAP14879 has found a total of 688 pores. 333 pores available for immediate sequencing April 16, 16:34
- Performing Mux Scan April 16, 16:32
- Mux scan for flow cell FAP14879 has found a total of 794 pores. 404 pores available for immediate sequencing April 16, 15:01
- Performing Mux Scan April 16, 14:59
- Mux scan for flow cell FAP14879 has found a total of 806 pores. 394 pores available for immediate sequencing April 16, 13:28
- Performing Mux Scan April 16, 13:25
- Mux scan for flow cell FAP14879 has found a total of 856 pores. 422 pores available for immediate sequencing April 16, 11:54
- Performing Mux Scan April 16, 11:52
- Mux scan for flow cell FAP14879 has found a total of 835 pores. 386 pores available for immediate sequencing April 16, 10:21
- Performing Mux Scan April 16, 10:18
- Mux scan for flow cell FAP14879 has found a total of 854 pores. 391 pores available for immediate sequencing April 16, 08:47
- Performing Mux Scan April 16, 08:45
- Mux scan for flow cell FAP14879 has found a total of 888 pores. 410 pores available for

- immediate sequencing April 16, 07:14
- Performing Mux Scan April 16, 07:11
- Mux scan for flow cell FAP14879 has found a total of 910 pores. 414 pores available for immediate sequencing April 16, 05:40
- Performing Mux Scan April 16, 05:38
- Mux scan for flow cell FAP14879 has found a total of 976 pores. 460 pores available for immediate sequencing April 16, 04:07
- Performing Mux Scan April 16, 04:04
- Mux scan for flow cell FAP14879 has found a total of 980 pores. 447 pores available for immediate sequencing April 16, 02:33
- Performing Mux Scan April 16, 02:31
- Mux scan for flow cell FAP14879 has found a total of 1013 pores. 470 pores available for immediate sequencing April 16, 01:00
- Performing Mux Scan April 16, 00:57
- Mux scan for flow cell FAP14879 has found a total of 1019 pores. 467 pores available for immediate sequencing April 15, 23:26
- Performing Mux Scan April 15, 23:24
- Mux scan for flow cell FAP14879 has found a total of 1021 pores. 470 pores available for immediate sequencing April 15, 21:53
- Performing Mux Scan April 15, 21:51
- Mux scan for flow cell FAP14879 has found a total of 1054 pores. 476 pores available for immediate sequencing April 15, 20:19
- Performing Mux Scan April 15, 20:17
- Mux scan for flow cell FAP14879 has found a total of 1050 pores. 470 pores available for immediate sequencing April 15, 18:46
- Performing Mux Scan April 15, 18:44
- Mux scan for flow cell FAP14879 has found a total of 1051 pores. 473 pores available for immediate sequencing April 15, 17:13
- Performing Mux Scan April 15, 17:10
- Mux scan for flow cell FAP14879 has found a total of 1060 pores. 476 pores available for immediate sequencing April 15, 15:39
- Performing Mux Scan April 15, 15:37
- Mux scan for flow cell FAP14879 has found a total of 1078 pores. 476 pores available for immediate sequencing April 15, 14:06
- Performing Mux Scan April 15, 14:03
- Mux scan for flow cell FAP14879 has found a total of 1081 pores. 477 pores available for immediate sequencing April 15, 12:32
- Performing Mux Scan April 15, 12:30
- Mux scan for flow cell FAP14879 has found a total of 1100 pores. 487 pores available for immediate sequencing April 15, 10:59
- Performing Mux Scan April 15, 10:56
- Starting sequencing procedure April 15, 10:56
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C April 15, 10:53