



## Run Info

Host Name	GXB01190 (localhost)
Experiment Name	EIMock_2kbp_NoEnrichRAD_050121
Sample ID	EIMock_2kbp_NoEnrichRAD_050121
Run ID	7f4a7183-0cc2-4ba5-847e-a067f8b22fa5
Flow Cell Id	FAO53362
Start Time	January 5, 12:56
Run Length	56m

## Run Summary

Reads Generated	229.07 K
Passed Bases	361.42 Mb
Failed Bases	23.1 Mb
Estimated Bases	414.79 Mb


## Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-RAD004
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=["/data/references/the7references.fasta"]
Read Filtering	min_qscore=7


## Versions

MinKNOW	20.10.6
MinKNOW Core	4.1.2
Bream	6.1.4
Guppy	4.2.3

Cumulative Output Reads




Cumulative Output Bases




**Read Length Histogram Estimated Bases - Outliers Discarded**

Estimated N50: 2.44 K




**Read Length Histogram Basecalled Bases - Outliers Discarded**

Estimated N50: 2.42 K




**Read Length Histogram Estimated Bases**

Estimated N50: 2.44 K




**Read Length Histogram Basecalled Bases**


Estimated N50: 2.42 K




Alignment Target Hits (reads)




Alignment Target Hits (bases)




Duty Time Grouped




Duty time Categorised




### Mux Scan Grouped




### Mux Scan Categorised



**Temperature History.**




**Bias Voltage History.**







Translocation Speed



QScore



**Disk Write Performance**



### **Run Debug Messages**

- Mux scan for flow cell FAO53362 has found a total of 1048 pores. 475 pores available for immediate sequencing January 5, 13:02
- Performing Mux Scan January 5, 13:00
- Starting sequencing procedure January 5, 13:00
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C January 5, 12:56