



Run Info

Host Name	GXB03020 (localhost)
Experiment Name	ReadUntil_38kbp_LowtoHigh_con_15042021
Sample ID	ReadUntil_38kbp_LowtoHigh_con_15042021
Run ID	9f8e23e9-1c06-43c2-9a43-5b85998e7fc2
Flow Cell Id	FAP14753
Start Time	April 15, 10:25
Run Length	1h 0m

Run Summary

Reads Generated	45.63 K
Passed Bases	412.17 Mb
Failed Bases	25.69 Mb
Estimated Bases	432.72 Mb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=["/data/the7references.fasta"]
Read Filtering	min_qscore=7


Versions

MinKNOW	21.02.5
MinKNOW Core	4.2.5
Bream	6.1.10
Guppy	4.3.4

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 35.47 K




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 35.02 K




Read Length Histogram Estimated Bases

Estimated N50: 35.47 K




Read Length Histogram Basecalled Bases


Estimated N50: 35.02 K




Alignment Target Hits (reads)




Alignment Target Hits (bases)




Duty Time Grouped




Duty time Categorized




Mux Scan Grouped




Mux Scan Categorised




Temperature History.




Bias Voltage History.




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAP14753 has found a total of 1074 pores. 467 pores available for immediate sequencing April 15, 10:31
- Performing Mux Scan April 15, 10:28
- Starting sequencing procedure April 15, 10:28
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C April 15, 10:25