



Run Info

Host Name	GXB01190 (localhost)
Experiment Name	ReadUntilMock_15kbSE_NoEnrich_08122020
Sample ID	ReadUntilMock_15kbSE_NoEnrich_08122020
Run ID	194daf26-2b0b-4784-a335-6548f883d3fb
Flow Cell Id	FAO52906
Start Time	December 8, 11:01
Run Length	38m

Run Summary

Reads Generated	18.4 K
Passed Bases	218.15 Mb
Failed Bases	7.85 Mb
Estimated Bases	230.36 Mb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Enabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Alignment	reference_files=["/data/references/the7references.fasta"]
Read Filtering	min_qscore=7


Versions

MinKNOW	20.10.6
MinKNOW Core	4.1.2
Bream	6.1.4
Guppy	4.2.3

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 26.14 K




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 25.85 K




Read Length Histogram Estimated Bases

Estimated N50: 26.14 K




Read Length Histogram Basecalled Bases


Estimated N50: 25.85 K




Alignment Target Hits (reads)




Alignment Target Hits (bases)




Duty Time Grouped




Duty time Categorized




Mux Scan Grouped




Mux Scan Categorised




Temperature History.




Bias Voltage History.




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAO52906 has found a total of 1204 pores. 496 pores available for immediate sequencing December 8, 11:07
- Performing Mux Scan December 8, 11:05
- Starting sequencing procedure December 8, 11:05
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C December 8, 11:01